

Product Datasheet

Rat SPD ELISA Kit (orb442329)

Description	Rat Surfactant Associated Protein D (SPD) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, lung lavage fluid and other biological fluids.
Sample Types	serum, plasma, lung lavage fluid and other biological fluids.
Uniprot ID	P35248
Sensitivity	34pg/mL
Entrez	25350
Expiration Date	Please enquire.