

Product Datasheet

Mouse SPINK5 ELISA Kit (orb442338)

Description	Mouse Serine Peptidase Inhibitor Kazal Type 5 (SPINK5) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q148R4
Sensitivity	28pg/mL
Entrez	72432
Expiration Date	Please enquire.