

Product Datasheet

Human SPINT2 ELISA Kit (orb442339)

Description	Human Serine Peptidase Inhibitor Kunitz Type 2 (SPINT2) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O43291
Sensitivity	0.103ng/mL
Entrez	10653
Expiration Date	Please enquire.