

Product Datasheet

Human SPLUNC2 ELISA Kit (orb442342)

Description	Human Short Palate, Lung And Nasal Epithelium Carcinoma Associated Protein 2 (SPLUNC2) Sandwich ELIS...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q96DR5
Sensitivity	0.059ng/mL
Entrez	140683
Expiration Date	Please enquire.