

Product Datasheet

Human SPLUNC3 ELISA Kit (orb442343)

Description	Human Short Palate, Lung And Nasal Epithelium Carcinoma Associated Protein 3 (SPLUNC3) Sandwich ELIS...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9BQP9
Sensitivity	0.066ng/mL
Entrez	128861
Expiration Date	Please enquire.