

Product Datasheet

Human ST14 ELISA Kit (orb442377)

Description	Human Suppression Of Tumorigenicity 14 (ST14) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9Y5Y6
Sensitivity	0.058ng/mL
Entrez	6768
Expiration Date	Please enquire.