

Product Datasheet

Mouse STAT3 ELISA Kit (orb442384)

Description	Mouse Signal Transducer And Activator Of Transcription 3 (STAT3) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P42227
Sensitivity	0.31ng/mL
Entrez	20848
Expiration Date	Please enquire.