

Product Datasheet

Human STAT5A ELISA Kit (orb442387)

Description	Human Signal Transducer And Activator Of Transcription 5A (STAT5A) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P42229
Sensitivity	0.112ng/mL
Entrez	6776
Expiration Date	Please enquire.