

Product Datasheet

Mouse STC1 ELISA Kit (orb442393)

Description	Mouse Stanniocalcin 1 (STC1) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O55183
Sensitivity	29.3pg/mL
Entrez	20855
Expiration Date	Please enquire.