

Product Datasheet

Human STEAP2 ELISA Kit (orb442396)

Description	Human Six Transmembrane Epithelial Antigen Of The Prostate 2 (STEAP2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q8NFT2
Sensitivity	0.062ng/mL
Entrez	261729
Expiration Date	Please enquire.