

Product Datasheet

FLNA (phospho-Tyr1046) antibody [Out of stock] (orb4424)

Description	Rabbit polyclonal to FLNA (phospho-Tyr1046)
Species/Host	Rabbit
Conjugation	Unconjugated
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human FLNA around the phosphorylation site of Tyr1046:VT(p-Y)DG
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	291kDa
Uniprot ID	Q8BTM8, P21333
Dilution Range	ELISA: 1:500-1000, IHC-P: 1:100-500, IF: 1:100-500
Expiration Date	12 months from date of receipt.