

Product Datasheet

Bovine SUOX ELISA Kit (orb442413)

Description	Bovine Sulfite Oxidase (SUOX) Sandwich ELISA Kit
Reactivity	Bovine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q3MHX0
Sensitivity	0.115ng/mL
Entrez	509837
Expiration Date	Please enquire.