

Product Datasheet

Human TAF13 ELISA Kit (orb442437)

Description	Human TATA Box Binding Protein Associated Factor 13 (TAF13) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q15543
Sensitivity	0.128ng/mL
Entrez	6884
Expiration Date	Please enquire.