

Product Datasheet

Rat TAFI ELISA Kit (orb442440)

Description	Rat Thrombin Activatable Fibrinolysis Inhibitor (TAFI) Sandwich ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	Q9EQV9
Sensitivity	0.56ng/mL
Entrez	113936
Expiration Date	Please enquire.