

Product Datasheet

Human TAP1 ELISA Kit (orb442446)

Description	Human Antigen Peptide Transporter 1 (TAP1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q03518
Sensitivity	0.134ng/mL
Entrez	6890
Expiration Date	Please enquire.