

Product Datasheet

Rat TAUT ELISA Kit (orb442459)

Description	Rat Taurine Transporter (TAUT) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P31643
Sensitivity	0.052ng/mL
Entrez	29464
Expiration Date	Please enquire.