

Product Datasheet

Human TBX21 ELISA Kit (orb442467)

Description	Human T-Box Protein 21 (TBX21) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	cell lysates and other biological fluids.
Sample Types	cell lysates and other biological fluids.
Uniprot ID	Q9UL17
Sensitivity	0.052ng/mL
Entrez	30009
Expiration Date	Please enquire.