

Product Datasheet

Human TECK ELISA Kit (orb442477)

Description	Human Thymus Expressed Chemokine (TECK) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O15444
Sensitivity	11.5pg/mL
Entrez	6370
Expiration Date	Please enquire.