


Product Datasheet

CD300a antibody (orb44248)

Description	Mouse Monoclonal to CD300a.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	HPB human acute lymphoid leukemia cell line
Target	CD300a
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2 µg/ml. Western blotting: Positive material: HPB cell line, negative material: JURKAT human cell line, RAJI human cell line; non-reducing conditions.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q9UGN4
Entrez	11314
Dilution Range	Flow cytometry: Recommended dilution: 2 µg/ml. Western blotting: Positive material: HPB cell line, negative material: JURKAT human cell line, RAJI human cell line; non-reducing conditions.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...