

Product Datasheet

Human TERF2IP ELISA Kit (orb442481)

Description	Human Telomeric Repeat Binding Factor 2 Interacting Protein (TERF2IP) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9NYB0
Sensitivity	0.66ng/mL
Entrez	54386
Expiration Date	Please enquire.