

Product Datasheet

Rat TFF1 ELISA Kit (orb442488)

Description	Rat Trefoil Factor 1 (TFF1) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q63467
Sensitivity	27pg/mL
Entrez	117270
Expiration Date	Please enquire.