

Product Datasheet

Human TFF3 ELISA Kit (orb442492)

Description	Human Trefoil Factor 3, Intestinal (TFF3) Sandwich ELISA Kit
Reactivity	Human
Range	125-8000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q07654
Sensitivity	46pg/mL
Entrez	7033
Expiration Date	Please enquire.