

Product Datasheet

Rat TGM2 ELISA Kit (orb442516)

Description	Rat Transglutaminase 2, Tissue (TGM2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10mU/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9WVJ6
Sensitivity	0.042mU/mL
Entrez	56083
Expiration Date	Please enquire.