

Product Datasheet

Human TGM4 ELISA Kit (orb442519)

Description	Human Transglutaminase 4, Prostate (TGM4) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P49221
Sensitivity	0.33ng/mL
Entrez	7047
Expiration Date	Please enquire.