

Product Datasheet

Bovine THBS1 ELISA Kit (orb442524)

Description	Bovine Thrombospondin 1 (THBS1) Sandwich ELISA Kit
Reactivity	Bovine
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum and plasma.
Sample Types	serum and plasma.
Uniprot ID	Q28178
Sensitivity	22pg/mL
Entrez	281530
Expiration Date	Please enquire.