

Product Datasheet

Mouse THBS2 ELISA Kit (orb442529)

Description	Mouse Thrombospondin 2 (THBS2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma, tissue homogenates and other biological fluids.
Sample Types	plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q03350
Sensitivity	0.124ng/mL
Entrez	21826
Expiration Date	Please enquire.