

Product Datasheet

Mouse THBS4 ELISA Kit (orb442533)

Description	Mouse Thrombospondin 4 (THBS4) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9Z1T2
Sensitivity	0.125ng/mL
Entrez	21828
Expiration Date	Please enquire.