

Product Datasheet

Human TIMD4 ELISA Kit (orb442556)

Description	Human T-Cell Immunoglobulin And Mucin Domain Containing Protein 4 (TIMD4) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q96H15
Sensitivity	0.062ng/mL
Entrez	91937
Expiration Date	Please enquire.