

Product Datasheet

Porcine TIMP1 ELISA Kit (orb442562)

Description	Porcine Tissue Inhibitors Of Metalloproteinase 1 (TIMP1) Sandwich ELISA Kit
Reactivity	Porcine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P35624
Sensitivity	13.6pg/mL
Entrez	396862
Expiration Date	Please enquire.