

Product Datasheet

Canine TIMP2 ELISA Kit (orb442566)

Description	Canine Tissue Inhibitors Of Metalloproteinase 2 (TIMP2) Sandwich ELISA Kit
Reactivity	Canine
Range	6.25-400ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9TTY1
Sensitivity	1.5ng/mL
Entrez	403633
Expiration Date	Please enquire.