

Product Datasheet

Equine TIMP2 ELISA Kit (orb442567)

Description	Equine Tissue Inhibitors Of Metalloproteinase 2 (TIMP2) Sandwich ELISA Kit
Reactivity	Equine
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	077717
Sensitivity	1.49ng/mL
Entrez	100034134
Expiration Date	Please enquire.