

Product Datasheet

Chicken TIMP3 ELISA Kit (orb442574)

Description	Chicken Tissue Inhibitors Of Metalloproteinase 3 (TIMP3) Sandwich ELISA Kit
Reactivity	Gallus
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P26652
Sensitivity	0.076ng/mL
Entrez	396483
Expiration Date	Please enquire.