

Product Datasheet

Porcine TIMP4 ELISA Kit (orb442587)

Description	Porcine Tissue Inhibitors Of Metalloproteinase 4 (TIMP4) Sandwich ELISA Kit
Reactivity	Porcine
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9TTB6
Sensitivity	0.068ng/mL
Entrez	396774
Expiration Date	Please enquire.