

Product Datasheet

Rat TIMP4 ELISA Kit (orb442588)

Description	Rat Tissue Inhibitors Of Metalloproteinase 4 (TIMP4) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P81556
Sensitivity	24pg/mL
Entrez	680130
Expiration Date	Please enquire.