

Product Datasheet

Rabbit TIMP4 ELISA Kit (orb442589)

Description	Rabbit Tissue Inhibitors Of Metalloproteinase 4 (TIMP4) Sandwich ELISA Kit
Reactivity	Rabbit
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	G1TAK3
Sensitivity	36pg/mL
Entrez	100008691
Expiration Date	Please enquire.