


Product Datasheet

CD30 antibody (orb44260)

Description	Mouse Monoclonal to CD30.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	Expression vector containing CD30 cDNA (booster suspension of THP-1 cell line)
Target	CD30
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 5 µg/ml. Immunohistochemistry (paraffin sections): Recommended dilution: 5 µg/ml, positive tissue: thymus.
Isotype	Mouse IgG2a kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P28908
Entrez	943
Dilution Range	Flow cytometry: Recommended dilution: 5 µg/ml. Immunohistochemistry (paraffin sections): Recommended dilution: 5 µg/ml, positive tissue: thymus.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...