

Product Datasheet

Mouse TLR3 ELISA Kit (orb442610)

Description	Mouse Toll Like Receptor 3 (TLR3) Sandwich ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q99MB1
Sensitivity	23.64pg/mL
Entrez	142980
Expiration Date	Please enquire.