

Product Datasheet

Human TLR8 ELISA Kit (orb442620)

Description	Human Toll Like Receptor 8 (TLR8) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9NR97
Sensitivity	0.108ng/mL
Entrez	51311
Expiration Date	Please enquire.