

Product Datasheet

Canine Thrombomodulin ELISA Kit (orb442622)

Description	Canine Thrombomodulin (TM) Sandwich ELISA Kit
Reactivity	Canine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	Q5W7P8
Sensitivity	5.8pg/mL
Entrez	474355
Expiration Date	Please enquire.