

Product Datasheet

Rat TMEM27 ELISA Kit (orb442625)

Description	Rat Transmembrane Protein 27 (TMEM27) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9ESG3
Sensitivity	0.046ng/mL
Entrez	57395
Expiration Date	Please enquire.