

Product Datasheet

Rat Thrombomodulin ELISA Kit (orb442635)

Description	Rat Thrombomodulin (TM) Sandwich ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	O35370
Sensitivity	0.57ng/mL
Entrez	83580
Expiration Date	Please enquire.