

Product Datasheet

Mouse TNNI3 ELISA Kit (orb442650)

Description	Mouse Troponin I Type 3, Cardiac (TNNI3) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P48787
Sensitivity	12.7pg/mL
Entrez	21954
Expiration Date	Please enquire.