

Product Datasheet

Rabbit TNNT2 ELISA Kit (orb442660)

Description	Rabbit Troponin T Type 2, Cardiac (TNNT2) Sandwich ELISA Kit
Reactivity	Rabbit
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P09741
Sensitivity	6.2pg/mL
Entrez	100009428
Expiration Date	Please enquire.