

Product Datasheet

Human TOMM70A ELISA Kit (orb442668)

Description	Human Translocase Of Outer Mitochondrial Membrane 70A (TOMM70A) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9CZW5
Sensitivity	0.101ng/mL
Entrez	28185
Expiration Date	Please enquire.