

Product Datasheet

Human TP63 ELISA Kit (orb442680)

Description	Human Tumor Protein P63 (TP63) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9H3D4
Sensitivity	0.057ng/mL
Entrez	8626
Expiration Date	Please enquire.