

Product Datasheet

Rat TPH1 ELISA Kit (orb442689)

Description	Rat Tryptophan Hydroxylase 1 (TPH1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P09810
Sensitivity	0.116ng/mL
Entrez	24848
Expiration Date	Please enquire.