

Product Datasheet

Mouse Thymidine Phosphorylase ELISA Kit (orb442700)

Description	Mouse Thymidine Phosphorylase (TP) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q99N42
Sensitivity	0.142ng/mL
Entrez	72962
Expiration Date	Please enquire.