

Product Datasheet

Bovine TPO ELISA Kit (orb442701)

Description	Bovine Thrombopoietin (TPO) Sandwich ELISA Kit
Reactivity	Bovine
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	C8CJH1
Sensitivity	31pg/mL
Entrez	525013
Expiration Date	Please enquire.