

Product Datasheet

Mouse TPP1 ELISA Kit (orb442707)

Description	Mouse Tripeptidyl Peptidase I (TPP1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O89023
Sensitivity	0.055ng/mL
Entrez	12751
Expiration Date	Please enquire.