

Product Datasheet

Human TRAF3 ELISA Kit (orb442722)

Description	Human TNF Receptor Associated Factor 3 (TRAF3) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q13114
Sensitivity	29pg/mL
Entrez	7187
Expiration Date	Please enquire.